

Protect Cleaning Indicators
Consistent Monitoring



Product Description

The **Protect Cleaning Indicator (PCI)** is a verification test used to monitor the effectiveness of the cleaning process in washer / disinfectors.

The indicator mimics soiled surgical instruments, providing an appropriate challenge for all factors in automated cleaning including time, temperature, pressure and detergent, including enzyme activity.

Indicator Coupon for Cleaning Performance Test

Designed for routine control of cleaning processes and for monitoring automatic washers.

Holder

Designed to mimic the cleaning challenges posed by a large process load during the washing cycle.

Product Information

Description	Product Number	Dimensions	Quantity #
Indicator Coupon	1108326	1.25 x 3 in	200 per package
Holder (Blue)	1108327	1.25 x 1 x .25 in	1 per package

Product Features

- The Protect Cleaning Indicator was thoroughly tested to ensure the coupon and holder provide a sufficient challenge and consistent monitoring results.
- Clips to the side of baskets or trays during testing, providing a better challenge for blind spots in the washer / disinfectors.
- PCI is made of a synthetic recyclable substrate which is more sustainable and cost-effective than metal-based indicators.
- PCI has a three-year shelf life, more than double the shelf life of other cleaning verification products.